## SAFE® ACCESS TUBE REFORMAT 96

## **Automated Tube Carrier**

- Fully automated tube carrier
- Designed for the bulk transfer of sample tubes from a filled SBS/ANSI 96 tube rack to an empty rack
- TUBE REFORMAT 96 works with most tube racks
- The contents of the rack (containing 1 to 96 tubes) can be transferred to the same locations in the receiving rack at a rate of approximately 1 transfer per minute
- Easy to operate
- Small footprint
- Dimensions: 245 mm x 295 mm x 388 mm

Item Number

RF96

Order Information



