



CERTIFICATE OF ANALYSIS

ESR_COA_SI305#2088B.pdf

| | | | |
|-------------------------------------|--|-------------------------------------|---------------------------------|
| Product Code SI 305.300-A | Name Latex Controls (30 Tests) | Quantity manufactured 178 | Type of Tubes Greiner |
|-------------------------------------|--|-------------------------------------|---------------------------------|

| | | | | |
|----------------------------|---|----------------------------------|--|--|
| Lot Number 2088B | Manufacturing Date 05/05/2020 | Expiry Date 01/11/2020 | | |
|----------------------------|---|----------------------------------|--|--|

| Standard Values | | | Measured Values | | | | Internal Aggregation | Volume checked | | |
|-----------------|-----------------------------|---------|-----------------------------|-----------|----|--------------|----------------------|-------------------------------------|----|-----|
| Transmittance | ESR (mm/h) (min-Val-max) | | ESR (mm/h) (min-Val-max) | | | Instrument | Serial Number | OK | | |
| Level 2 | 1534 | 6 - 11 | 9 | 9 | 9 | TEST 1 BCL | 2173 | <input checked="" type="checkbox"/> | NO | YES |
| Level 3 | 980 | 15 - 22 | 19 | 19 | 20 | TEST1 THL | T13290 | <input checked="" type="checkbox"/> | NO | YES |
| Level 4 | 637 | 56 - 74 | 65 | 66 | 68 | ROLLER 20 PN | 133 | <input checked="" type="checkbox"/> | NO | YES |

Controlled By **Cirio Alberto** Date of Control **05/05/2020**
 Approved By Manufacturing Manager **Pividori Denis** Date of Approval **05/05/2020**

We certify that the lot in object has been manufactured and tested according with ALIFAX manufacturing and quality control procedures and specifications. These requirements have been established according to the European Directive 98/79/EC. We confirm that the above referenced product conformed to the specifications and then it is released.

Released By Quality **Viel Sara** Date of Release **06/05/2020**

This is a computer-generated document. No signature is required.

Registered Office

Via Petrarca, 2/1 - 35020 Polverara (PD) Italy
 Tel (+39) 049 0992000 - Fax (+39) 049 5855434
 Mail PEC: info.alifax@pec.alifax.com
info@alifax.com - www.alifax.com



Manufacturing Site

Via Merano, 30 - 33045 Nimis (UD) Italy
 Tel (+39) 0432 547454 - Fax (+39) 0432 547378
 Mail PEC: info.udine@pec.alifax.com
info.udine@alifax.com