

Product Code	Name	Quantity manufactured	
<b>SI 390.900</b>	<b>URO-QUICK SCREENING KIT 360 TEST</b>	<b>26</b>	

Lot Number	Manufacturing Date	Expiry Date	CAP color	Printed Band
<b>SCR14179</b>	<b>24/07/2014</b>	<b>24/09/2015</b>	<b>Green</b>	<b>White</b>

## Eugonic Broth

Lot Number **SCR14179**  
 Expiry Date **24/09/2015**  
 CAP color **Green**  
 Printed Band **White**

	Physical parameters <i>(Volume, transparency, colour)</i>	Decontamination at 37°C after 72 hours	Growth inter batch	Growth intra batch
Expected	Compliant	No growth	CV% <= 35%	CV% <= 35%
Found	OK	OK	7.7%	< 5%

Check of the presence of the magnetic stirrer in each vial contained in the kit: **YES**

Storage Conditions: **min 4°C - Max 30°C**

Compliance: **YES**

Date of Control: **31/07/2014**

Date of Approval: **01/08/2014**

Tested: **Cirio Alberto**

Quality Control Manager: **Cossettini Marco**


