

CERTIFICATE OF ANALYSIS

MIC_COA_SI390900#SCR14180.pdf

Product Code SI 390.900	Name URO-QUICK SCREENING KIT 360 TEST	Quantity manufactured 26	
Lot Number SCR14180	Manufacturing Date 25/07/2014	Expiry Date 25/09/2015	CAP color Green Printed Band White

Eugonic Broth

Lot Number **SCR14180**
 Expiry Date **25/09/2015**
 CAP color **Green**
 Printed Band **White**

	Physical parameters <i>(Volume, transparency, colour)</i>	Decontamination at 37°C after 72 hours	Growth inter batch	Growth intra batch
Expected	Compliant	No growth	CV% ≤ 35%	CV% ≤ 35%
Found	OK	OK	12.5%	17.7%

Check of the presence of the magnetic stirrer in each vial contained in the kit: **YES**

Storage Conditions: **min 4°C - Max 30°C**

Compliance: **YES**

Date of Control **31/07/2014**

Date of Approval **01/08/2014**

Tested **Cirio Alberto**

Quality Control Manager **Cossettini Marco**


