

---

### Certificate of Quality Control

Manufacturer	<b>SIFIN</b> Institut für Immunpräparate und Nährmedien GmbH Berlin, Berliner Allee 317-321, 13088 Berlin, Germany
Product	<b>Enteroclon Anti-Salmonella O 5</b>
Catalogue No. (Package size)	TR 1303 (1 ml)
Lot No.	<b>150 11 11</b>
Expiry date	<b>2016-11</b>
Intended use	Enteroclon Anti-Salmonella O 5 records all Salmonella strains of group B, as long as they contain antigen O 5.
Specificity	No agglutination with heterologous Salmonella O antigens and Salmonella H antigens.
Preservatives	Sodium azide 0.9 mg / ml
Storage	2...8 °C

Berlin, 17-02-12

  
Head of Quality Control