


Certificate of Quality Control

Manufacturer	SIFIN Institut für Immunpräparate und Nährmedien GmbH Berlin, Berliner Allee 317-321, 13088 Berlin, Germany
Product	Enteroclon Anti-Salmonella O 41
Catalogue No. (Package size)	TR 1347 (1 ml)
Lot No.	161 12 11
Expiry date	2016-12
Intended use	Enteroclon Anti-Salmonella O 41 records all Salmonella strains, as long as they contain antigen O 41.
Specificity	No other agglutinations with heterologous Salmonella O antigens and Salmonella H antigens are known.
Preservatives	Sodium azide 0.9 mg / ml
Storage	2...8 °C

Berlin,

17.04.12


Head of Quality Control