

Certificate of Quality Control

Manufacturer	SIFIN Institut für Immunpräparate und Nährmedien GmbH Berlin, Berliner Allee 317-321, 13088 Berlin, Germany
Product	Enteroclon Anti-Salmonella H g,m
Catalogue No. (Package size)	TR 1408 (1 ml)
Lot No.	148 09 13
Expiry date	2015-09
Intended use	Enteroclon Anti-Salmonella H g,m records all Salmonella strains, as long as they contain in an antigen combination the H g and the H m antigens (e.g. f,g; f,g,s; f,g,t; g,s,t; g,q; g,t; g,z ₅₁ or g,m; g,m,s,t; m,t) in slide agglutination.
Specificity	No agglutination with heterologous Salmonella H antigens and Salmonella O antigens.
Preservatives	Sodium azide 0.9 mg/ ml
Storage	2...8 °C

Berlin, 28.10.13



Head of Quality Control