

## Certificate of Quality Control

CERTIFIED QUALITY MANAGEMENT SYSTEM IN  
ACCORDANCE WITH DIN EN ISO 9001 AND DIN EN ISO 13485

sifin diagnostics gmbh

Berliner Allee 317-321  
13088 Berlin, Germany

T: +49 30 927 030-0  
F: +49 30 927 030-30  
info@sifin.de  
www.sifin.de

<b>Product</b>	<b>Anti-Salmonella H:L</b>
Catalogue no.	TR 1412
Package size	1 ml
<b>Lot no.</b>	<b>78 04 14</b>
<b>Expiry date</b>	<b>2017-04</b>
Intended use	Anti-Salmonella H:L records all <i>Salmonella</i> strains, as long as they contain antigen H:l (e.g. H:l,v; H:l,w; H:l,z <sub>13</sub> ; H:l,z <sub>28</sub> ; H:l,z <sub>13</sub> ,z <sub>28</sub> ) in slide agglutination.
Specificity	No agglutination with heterologous <i>Salmonella</i> H antigens and <i>Salmonella</i> O antigens.
Preservatives	Sodium azide 0.9 mg/ml
Storage	2 °C... 8 °C

Berlin,  
place, date

10.06.14

signature Head of Quality Control