

## Certificate of Quality Control

CERTIFIED QUALITY MANAGEMENT SYSTEM IN  
ACCORDANCE WITH DIN EN ISO 9001 AND DIN EN ISO 13485

sifin diagnostics gmbh

Berliner Allee 317-321  
13088 Berlin, Germany

T: +49 30 927 030-0  
F: +49 30 927 030-30  
info@sifin.de  
www.sifin.de

<b>Product</b>	<b>Anti-Salmonella H:L</b>
Catalogue no.	TR 1412
Package size	1 ml
<b>Lot no.</b>	<b>255 12 15</b>
<b>Expiry date</b>	<b>2018-12</b>
Intended use	Anti-Salmonella H:L records all <i>Salmonella</i> strains, as long as they contain antigen H:I (e.g. H:I,v; H:I,w; H:I,Z <sub>13</sub> ; H:I,Z <sub>28</sub> ; H:I,Z <sub>13</sub> ,Z <sub>28</sub> ) in slide agglutination.
Specificity	No agglutination with heterologous <i>Salmonella</i> H antigens and <i>Salmonella</i> O antigens.
Preservatives	Sodium azide 0.9 mg/ml
Storage	2 °C...8 °C

Berlin,

04.02.16

place, date



signature Head of Quality Control