

Certificate of Quality Control

Manufacturer **SIFIN** Institut für Immunpräparate und Nährmedien GmbH Berlin,
Berliner Allee 317-321, 13088 Berlin, Germany

Product **Enteroclon Anti-Salmonella H y**

Catalogue No. TR 1423 (1 ml)
(Package size)

Lot No. **159 12 11**

Expiry date **2014-12**

Intended use Enteroclon Anti-Salmonella H y records all Salmonella strains, as long as they contain antigen H y.

Specificity No agglutination with heterologous Salmonella O antigens and Salmonella H antigens.

Preservatives Sodium azide 0.9 mg/ ml

Storage 2...8 °C

Berlin, 22.12.11


Head of Quality Control