

Certificate of Quality Control

CERTIFIED QUALITY MANAGEMENT SYSTEM IN ACCORDANCE WITH DIN EN ISO 9001 AND DIN EN ISO 13485

sifin diagnostics gmbh

Berliner Allee 317-321 13088 Berlin, Germany

T: +49 30 927 030-0 F: +49 30 927 030-30 info@sifin.de www.sifin.de

Product

Anti-Salmonella H:z₁₃

Catalogue no.

TR 1439

Package size

1 ml

Lot no.

66 03 14

Expiry date

2017-03

Intended use

Test Reagent Anti-Salmonella H:z₁₃ records all Salmonella strains, as long as

they contain antigen H:z₁₃ of the L complex in slide agglutination.

The product is not intended for use in vitro for the examination of samples

derived from the human body.

Specificity

No agglutination with Salmonella O antigens and heterologous Salmonella H

antigens, with the exception of H:1, H:z28 and H:z39.

Preservatives

Sodium azide 0,9 mg/ml

Storage

2 °C...8 °C

Berlin, 03,04,14

place, date

signature Head of Quality Control