

## Certificate of Quality Control

sifin diagnostics gmbh

CERTIFIED QUALITY MANAGEMENT SYSTEM IN  
ACCORDANCE WITH DIN EN ISO 9001 AND DIN EN ISO 13485

Berliner Allee 317-321  
13088 Berlin, Germany

T: +49 30 927 030-0  
F: +49 30 927 030-30  
info@sifin.de  
www.sifin.de

<b>Product</b>	<b>Anti-Shigella I</b> S. flexneri type 1-6, Gr.3,4 (y), 6, 7,8 (x), S. sonnei
Catalogue no.	TR 1811
Package size	1 ml
<b>Lot no.</b>	<b>1590217</b>
<b>Expiry date</b>	<b>2022-02</b>
Intended use	Anti-Shigella flexneri records all type- and group-specific antigens of <i>Shigella flexneri</i> in slide agglutination. as well as both forms of <i>Shigella sonnei</i> (Phase I-and Phase II).
Specificity	No agglutination with the serotypes of <i>Shigella dysenteriae</i> (type 1-10) and <i>Shigella boydii</i> (type 1 to 15).
Preservatives	Sodium azide 0.9 mg/ml
Storage	2 °C...8 °C

Berlin,  
place, date

25.09.2017

  
signature Head of Quality Control