

Certificate of Quality Control

CERTIFIED QUALITY MANAGEMENT SYSTEM IN
ACCORDANCE WITH DIN EN ISO 9001 AND DIN EN ISO 13485

sifin diagnostics gmbh

Berliner Allee 317-321
13088 Berlin, Germany

T: +49 30 700 144-0
F: +49 30 700 144-30
info@sifin.de
www.sifin.de

Product	Anti-Shigella I S. flexneri type 1-6, Gr.3,4 (y), 6, 7,8 (x), S. sonnei
Catalogue no.	TR1811
Package size	1 ml
Lot no.	5810720
Expiry date	2025-06
Intended use	Anti-Shigella flexneri records all type- and group-specific antigens of <i>Shigella flexneri</i> in slide agglutination. as well as both forms of <i>Shigella sonnei</i> (Phase I-and Phase II).
Specificity	No agglutination with the serotypes of <i>Shigella dysenteriae</i> (type 1-10) and <i>Shigella boydii</i> (type 1 to 15).
Preservatives	Sodium azide 0.9 mg/ml
Storage	2 °C...8 °C

Berlin,

14.09.2020

place, date

F. Ludwig

signature Head of Quality Control