


Certificate of Quality Control

Manufacturer	SIFIN Institut für Immunpräparate und Nährmedien GmbH Berlin, Berliner Allee 317-321, 13088 Berlin, Germany
Product	Test Reagent Anti-Coli III
Catalogue No. (Package size)	TS 2131 (1 ml)
Lot No.	06 01 12
Expiry date	2017-01
Intended use	Test Reagent Anti-Coli III records all E. coli strains, as long as they contain the antigens O25:(K11), O78:(K80), O103:(K-), O118:(K-), O124:(K72), O145:(K-), O157:(K-) and O164:(K-).
Specificity	No agglutination with E. coli strains which are recorded by Test Serum Anti-Coli I and Test Reagent Anti-Coli II.
Preservatives	Sodium azide 0.9 mg/ ml
Storage	2...8 °C

Berlin, 18.05.12


Head of Quality Control