


---

### Certificate of Quality Control

Manufacturer	<b>SIFIN</b> Institut für Immunpräparate und Nährmedien GmbH Berlin, Berliner Allee 317-321, 13088 Berlin, Germany
Product	<b>Test Reagent Anti-Coli III</b>
Catalogue No. (Package size)	TR 2131 (1 ml)
Lot No.	<b>68 07 10</b>
Expiry date	<b>2015-07</b>
Intended use	Test Reagent Anti-Coli III records all E. coli strains, as long as they contain the antigens O25:(K11), O78:(K80), O103:(K-), O118:(K-), O124:(K72), O145:(K-), O157:(K-) and O164:(K-).
Specificity	No agglutination with E. coli strains which are recorded by Test Serum Anti-Coli I and Test Reagent Anti-Coli II.
Preservatives	Sodium azide 0.9 mg/ ml
Storage	2...8 °C

Berlin, 17.12.10

  
\_\_\_\_\_  
Head of Quality Control