

Certificate of Quality Control

CERTIFIED QUALITY MANAGEMENT SYSTEM IN
ACCORDANCE WITH DIN EN ISO 9001 AND DIN EN ISO 13485

sifin diagnostics gmbh

Berliner Allee 317-321
13088 Berlin, Germany

T: +49 30 927 030-0
F: +49 30 927 030-30
info@sifin.de
www.sifin.de

Product	Anti-Coli III
Catalogue no.	TR 2131
Package size	1 ml
Lot no.	4990518
Expiry date	2023-05
Intended use	Anti-Coli III records all <i>E. coli</i> strains, as long as they contain the antigens O 25:(K 11), O 78:(K 80), O 103:(K-), O 118:(K-), O 124:(K 72), O 145:(K-), O 157:(K-) and O 164:(K-).
Specificity	No agglutination with <i>E. coli</i> strains which are recorded by Anti-Coli I and Anti-Coli II.
Preservatives	Sodium azide 0.9 mg/ ml
Storage	2 °C...8 °C

Berlin, 30.07.2018
place, date

J. Lindner
signature Head of Quality Control