

## Certificate of Quality Control

CERTIFIED QUALITY MANAGEMENT SYSTEM IN  
ACCORDANCE WITH DIN EN ISO 9001 AND DIN EN ISO 13485

sifin diagnostics gmbh

Berliner Allee 317-321  
13088 Berlin, Germany

T: +49 30 700 144-0  
F: +49 30 700 144-30  
info@sifin.de  
www.sifin.de

<b>Product</b>	<b>Anti-Salmonella H:t</b>
Catalogue no.	TS1418
Package size	1 ml
<b>Lot no.</b>	<b>4290621</b>
<b>Expiry date</b>	<b>2026-05</b>
Intended use	Anti-Salmonella H:t records all <i>Salmonella</i> strains, as long as they contain in an antigen combination the H:t (e. g. m,t; g,t; g,s,t; f,g,t; g,m,s,t) in slide agglutination.
Specificity	No agglutination with heterologous <i>Salmonella</i> H antigens and <i>Salmonella</i> O antigens.
Preservatives	Sodium azide 0.9 mg/ml
Storage	2 °C...8 °C

Berlin,  
place, date

26.08.21

signature Head of Quality/Control

7. Linder