

Certificate of Quality Control

Manufacturer **SIFIN** Institut für Immunpräparate und Nährmedien GmbH Berlin,
Berliner Allee 317-321, 13088 Berlin, Germany

Product **Test Serum Anti-Salmonella H_{Z₄Z₂₃}**

Catalogue No. (Package size) TS 1425 (1 ml)

Lot No. **138 08 12**

Expiry date **2017-09**

Intended use Test Serum Anti-Salmonella H_{Z₄Z₂₃} records all Salmonella strains,
as long as they contain antigens H_{Z₄Z₂₃} / H_{Z₄Z₂₄} / H_{Z₄Z₃₂}.

Specificity No agglutination with Salmonella O antigens and heterologous
Salmonella H antigens.

Preservatives Sodium azide 0.9 mg/ ml

Storage 2...8 °C

Berlin, 05.10.12


Head of Quality Control