

## Certificate of Quality Control

Manufacturer	<b>SIFIN</b> Institut für Immunpräparate und Nährmedien GmbH Berlin, Berliner Allee 317-321, 13088 Berlin, Germany
Product	<b>Test Serum Anti-Salmonella H z<sub>6</sub></b>
Catalogue No. (Package size)	TS 1426 (1 ml)
Lot No.	<b>78 06 11</b>
Expiry date	<b>2016-06</b>
Intended use	Enteroclon Anti-Salmonella H z <sub>6</sub> records all Salmonella strains, as long as they contain antigen H z <sub>6</sub> .
Specificity	No agglutination with heterologous Salmonella H antigens and Salmonella O antigens.
Preservatives	Sodium azide 0.9 mg/ ml
Storage	2...8 °C

Berlin, 25.07.11

R. R. K.  
Head of Quality Control