

• Certificate of Analysis

Product Name : MALDI Sepsityper Kit 50

Part No. 8270170

LOT No.: 0000290412

Use by¹⁾: 2019-01-31

Date of Production: 2018-02-12

Storage Conditions : All kit components must be stored at room temperature (20–25 °C).

<u>Test:</u>	<u>Specification</u>	<u>LOT Results</u>
Identification of <i>Candida glabrata</i> from a blood culture using the MALDI Biotyper and the MALDI Sepsityper Kit 50 following the kit sample preparation procedure.	Classification Result Report lists „ <i>Candida glabrata</i> ” as the best match. log(score) Value ≥ 2.0 of the best match.	Passed Passed
Identification of blank test samples based on one millilitre washing buffer using the MALDI Biotyper and the MALDI Sepsityper Kit 50 sample preparation procedure.	Classification Result Report lists „no peaks found” or “no reliable identification” as the best match. log(score) Value < 1.6 of the best match.	Passed

Remarks: ¹⁾Kit components should be handled carefully to avoid contaminations.

Bruker Daltonik GmbH, Germany
Fahrenheitstr. 4 – 28359 Bremen
phone +49-421-2205-0


 Bruker Daltonik GmbH
 Fahrenheitstraße 4
 D-28359 Bremen
 Tel: +49 (0) 421 - 2205-0
 Dr. Markus Kostrzewa, SAFETY OFFICER

Not for use in diagnostic procedures. For research use only.